

Disposable Gloves

The Fisherbrand™ range of single-use gloves are designed to provide maximum protection and comfort. They are the ideal gloves for everyday protection in the lab with new packaging designed to reduce waste and environmental impact.

Nitrile disposable gloves, comfort, PPE Category III

CE 0493

Choose a glove with enhanced comfort and tactile sensitivity. Fisherbrand comfort nitrile gloves are powder-free and latex-free with a beaded fit and textured fingers.

- Powder free, latex free
- Exceptional fit, feel and durability
- Ambidextrous with textured fingers and beaded cuff
- Non-sterile
- 200 gloves per standard box (size XL 180 gloves per box) reduces waste and environmental impact
- Thickness: 0.07mm (palm), 0.09mm (finger)
- Colour: blue
- Length 241.3mm
- Conforms to EN ISO 374/Type B



Ordering Information

Cat. No.	Size	Pack qty
15622367	XS	200
15632367	S	200
15642367	M	200
15652367	L	200
15662367	XL	180

Nitrile disposable gloves, extended cuff, PPE Category III

CE 0493

Choose a glove with exceptional durability and added protection. The Fisherbrand extended cuff nitrile gloves are extra-thick and fully textured for a superior grip.

- Powder free, latex free
- Fully textured finish for a superior wet or dry grip
- Ambidextrous and non-sterile with beaded cuff
- Thickness: 0.15mm (palm), 0.18mm (finger)
- Colour: teal green
- Length 305mm
- Conforms to EN ISO 374/Type B



Ordering Information

Cat. No.	Size	Pack qty
15612367	XS	50
15662357	S	50
15672357	M	50
15682357	L	50
15692357	XL	50

Nitrile disposable gloves with Aloe vera, PPE Category III

These gloves feature a textured surface to provide ideal grip for dry and wet working conditions. They are green for easy identification, and have a soothing Aloe vera coating on the inner surface making them ideal for prolonged use.

- Premium nitrile gloves offering the comfort and feel of latex
- Comfortable, strong and durable
- Aloe vera coating on the inner surface to protect hands from drying, chapping and cracking
- 100 ambidextrous gloves per standard box reduces waste and environmental impact
- PPE complies with Category III, 1.5 AQL, EN ISO 374/Type B
- Thickness 0.06mm



Ordering Information

Cat. No.	Size	Pack qty
15193269	S	100
15103279	M	100
11772779	L	100
15113279	XL	100

Latex disposable gloves with Aloe vera, PPE Category III

Fisherbrand latex examination gloves combine user comfort with provision of a protective barrier, and also contain soothing Aloe vera to help prevent your skin drying and chapping.

- Latex examination gloves providing a protective barrier
- Contain soothing Aloe vera to help prevent skin from drying, cracking and chapping
- Textured to provide excellent grip
- Powder free
- EN ISO 374/Type B (AQL 1.5)
- Length 240mm



Ordering Information

Cat. No.	Size	Pack qty
15153269	S	100
15173269	M	100
15163269	L	100
15183269	XL	100

© 2017 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

12567_EN