

Liquid chromatography

Thermo Scientific Vanquish Core Pre-configured bundle

Dependable workhorse for routine analysis

ThermoScientific™ Vanquish™ Core HPLC Systems deliver reliable, high-performance HPLC for routine workflows, combining robust hardware with intuitive control to ensure consistent, reproducible results. All systems have Platinum rated ACT 2.0 Ecolabels.

Designed for operational efficiency and seamless integration, Vanquish Core helps laboratories maximize uptime and analytical confidence.



Environmental
Performance Factor

71 Certified
July 2025



My Green Lab®
ACT® Ecolabel

Thermo Scientific™ Vanquish™ Core Pre-configuration

Cat. No.	Type	Description
31906275	System Base Vanquish Core	Foundational, stackable unit for the Vanquish Core HPLC system. The structural, electrical and fluidic core of the system supporting all modules
	Quaternary Pump C	Max pressure 700 bar. Flow Range 0-10 mL/min. Able to mix up to four different solvents for flexible gradient formation.
	Split Autosampler C	4 racks of 54 vials capacity with a volume injection range of 0.01-100uL
	Column Compartment C	Still and forced air heating for a temperature range of 5-85C. Column capacity of up to 4 columns of 100 mm length. Particle size 3-5 µm
31906276	Standard Flow Cell UHPLC	Standard 10 mm flow cell for Diode Array Detector. Stainless steel with 120 bar Software providing integrated instrument control, data acquisition, processing, and compliant-ready reporting.
31906274	Diode Array Detector	Acquire UV-Vis spectra in multiple channels over a wide dynamic range. Analyte must contain a chromophore.
31906284	1-day onsite training	One-day on site training to provide hands-on system operation guidance, method optimization assistance, and best-practice consultation.

Choose 1 power cord (country specific)

4 pack EU: 31906277 UK: 31906278 Swiss: 31906279 DK: 31906280



Thermo Scientific™ Chromeleon™ CDS 7.4 Pre-configuration

Cat. No.	Type	Description
31906283	Chromeleon Enterprise Client	Software providing integrated instrument control, data acquisition, processing, and compliant-ready reporting.
31906282	Chromeleon Enterprise Instrument	
31906281	Chromeleon License	

Optional additions to Pre-configuration

Cat. No.	Type	Description
31908733	computer	Dell XE4 w/ 27 inch monitor

The Thermo Scientific™ Vanquish™ Core HPLC pre-configuration bundle features a modular design that allows laboratories to tailor pump, sampler, column compartment, and detector configurations to meet changing analytical requirements, making it a versatile, general-purpose platform for diverse workflows. Seamless integration with Thermo Scientific™ Chromeleon™ Chromatography Data System (CDS) provides comprehensive instrument control, automation, data processing, and compliance-ready capabilities to streamline laboratory operations. Installation and 1 year warranty are included with Pre-configuration bundle purchase.

Contact your Fisher Scientific representative for further details and to verify if the preconfigured bundle is available in your region.

Beyond instrumentation and software configurations, we support an optional complete end-to-end solution, including columns, consumables, and sample preparation products, ensuring optimized performance across the full analytical workflow.



Optional Thermo Scientific™ consumable starters

Cat. No.	Type	Description	Qty.
10384323	HPLC column	Hypersil GOLD™ C18 Selectivity, 3 µm 100 x 4.6 mm	1
10364603		Hypersil GOLD™ aQ C18 Polar Endcapped, 3 µm, 100 x 4.6 mm	1
10690256		Hypersil GOLD™ Phenyl 3 µm, 100 x 4.6 mm	1
17324143	Vial	SureSTART™ 2 mL clear glass screw top vials, Level 1	100 pk.
17304133	Closure	SureSTART™ 9 mm blue AVCS screw PP caps, white silicone/red PTFE septum, Level 1	100 pk.
17393803	Vial + Closure Kit	SureSTART™ EPA Certified screw clear 12 mL vials, white polyethylene foam/white PTFE septum, white caps kit, Level 2	100 pk.
10362725	Syringe	National Target 3 mL Luer-Lok all plastic disposable syringe	100 pk.
10085634	Filter	Target2™ Nylon syringe filter, 0.25 µm, 30 mm, no prefilter	100 pk.
10269554		Target2™ PVDF syringe filter, 0.45 µm, 17 mm	100 pk.
10353015		Target2™ PTFE (Hydrophobic), 0.45 µm, 17 mm	100 pk.
10098942	SPE Cartridge	HyperSep™ SPE C18 cartridge, 500 mg, 6 mL	30 pk.

Find more information about our consumables at
eu.fishersci.com/go/chroma-column-selectionguide

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk **Germany:** fishersci.de **Ireland:** fishersci.ie
Italy: fishersci.it **Finland:** fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl **Norway:** fishersci.no
Portugal: fishersci.pt **Spain:** fishersci.es **Sweden:** fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

General Laboratory Equipment – Not For Diagnostic Procedures. © 2026 Thermo Fisher Scientific Inc. All rights reserved. Thermo Scientific™ general lab equipment, consumables, and software are not intended for use in diagnostic purposes in accordance with our product documentation, manuals, and labels. They are designated for General Laboratory Use Only. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. My Green Lab and ACT are trademarks owned and controlled by My Green Lab, Corp. This information is not intended to encourage use of these products in any manner that might infringe on the intellectual property rights of others. Not all products are available in all countries. Please consult your local sales representatives for details. **THER-2604-JS-0326, 5.26 update**